



Small Mask

Medium Mask

Mouthpiece

Large Mask

## The Monaghan Advantage

- ❖ Cost effective
- ❖ Anti-static, can be used directly out of package
- ❖ Increased aerosol suspension time
- ❖ Optimized delivery system for consistent medication availability
- ❖ Easy to clean - dishwasher safe
- ❖ Compatible with most commonly prescribed MDIs
- ❖ Not made or manufactured with BPA, natural rubber latex, phthalates, lead, or PVC

## MONAGHAN Z STAT® ANTI-STATIC VALVED HOLDING CHAMBER (aVHC)

The AEROCHAMBER® brand of valved holding chamber (aVHC) has been trusted by hospitals for over 35 years. Now you have the same Z STAT® aVHC found in hospitals available to you at home. The optimally designed anti-static chamber ensures consistent medication availability.

This cost-effective choice has patient-centric design features and helps optimize respirable dose of inhaler medications. A full range of products is available, designed specifically for each patient type.

## THE FACTS

- **Cost effective without sacrificing aerosol quality**
- **Dishwasher safe, convenient way to clean chamber and promote good cleaning habits.**
- **Consumer Friendly**
- **Easy to read user instructions**
- **Visually appealing design**
- **Made in the USA, product is manufactured in the US using US and Canadian components.**

MORE TIME  
TO BREATHE  
EASY




518.561.7330 phone  
800.833.9653 toll free  
518.561.5088 fax

Monaghan Medical Corporation  
Plattsburgh, New York 12901  
[monaghanmed.com](http://monaghanmed.com)

## THE FACTS (continued)

- **FLOWSIGNAL® Whistle to promote proper technique when transitioning from a mask to a mouthpiece device.**
- **Flow dynamic valve and baffle system prevents patients from exhaling into the chamber, ensuring the remaining medication is available for the next breath.**
- **Anti-static chamber provides consistent aerosol delivery right out of the package.**
- **Universal backpiece fits commonly prescribed inhalers for easy insertion and removal.**



## CLEANING IS EASY

- **MONAGHAN Z STAT® aVHC (12 month, reusable single patient use)**
  - Designed for 12 months of use at home, clinic, or hospital environment.
  - MONAGHAN Z STAT® aVHC chambers are dishwasher safe – follow care and use instructions included with the chamber.

## ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
79710Z	Monaghan Z STAT® aVHC with Mouthpiece	case/10
79750Z	Monaghan Z STAT® aVHC with Mouthpiece	case/50
88710Z	Monaghan Z STAT® aVHC with Small Mask	case/10
78710Z	Monaghan Z STAT® aVHC with Medium Mask	case/10
80710Z	Monaghan Z STAT® aVHC with Large Mask	case/10



Contact your Monaghan Sales Territory Manager for more information or a demonstration.

