



Can add a filter for added protection

Available with or without manometer

AccuPAP™

The AccuPAP Positive Airway Pressure device is an innovative treatment option used to provide positive airway pressure for lung volume expansion therapy. Indicated for the prevention and treatment of atelectasis, AccuPAP may aid in the reduction of air trapping associated with asthma and/or COPD as well as optimize delivery of bronchodilators.

- Adjustable pressure valve minimizes the need for manometer
- Treatment can be delivered with mouthpiece or mask attachment
- Design allows for convenient delivery of aerosol medication during PAP treatment
- Pressure port with cap and ambient air inlet
- Available with or without manometer
- Can add a filter for added protection

Item #	Description	Qty
301-6000	AccuPAP Positive Airway Pressure Device, Mouthpiece, 7' Tubing	10/bx
301-6001	AccuPAP Positive Airway Pressure Device, Mouthpiece, 7' Tubing, with Manometer	10/bx
301-6002	AccuPAP Positive Airway Pressure Device with Universal O ₂ Connector	10/bx
301-6003	AccuPAP Positive Airway Pressure Device with Manometer and Universal O ₂ Connector	10/bx
301-6200	APEP PEP Therapy System	10/bx