



## Stoma-Gel™



The Curaplex<sup>®</sup> Stoma-Gel™ was created to enhance patient comfort and decrease soreness for those early tracheostomy patients that experience discomfort due to pressure.

The Stoma-Gel™ is ergonomically designed to cushion the area between the tracheostomy tube flange and stoma site.

The unique and patented design allows for placement without removing the tracheostomy tube, making it convenient for clinician and comfortable for the patient.

The Stoma-Gel™ also features a revolutionary contoured and “rib-like” design, replicating shock absorbers, which create a comfortable interface for the tracheostomy tube and patient.

They are single patient use and individually packaged for convenience.

- Hypoallergenic and made from a medical-grade, mineral oil-based gel
- Unique, patented “pressure absorbing ribs” increase comfort and decrease soreness
- One size fits most- slight cut can be made to fit size 9 tracheostomy tubes or larger
- Replacement after 3 days of use is recommended
- No adhesive necessary
- Latex-Free

item #	description	qty
301-0903	Stoma-Gel™ Skin Barrier	50/bx